



COVID-19 IgG/IgM Antibody Rapid Test



- Whole blood (Finger prick available), serum and plasma sample
- Easy to use - no analyser required
- Store at room temperature or refrigerator



A range of rapid tests for the qualitative detection of antibodies (IgG and IgM) against COVID-19 (SARS-CoV-2) in whole blood, plasma or serum.

Product Features :

- ✓ Tests can be performed using whole-blood
- ✓ Visual results available in only 10 minutes
- ✓ 20 complete tests in one kit
- ✓ 96% accuracy - British-made reliability



The COVID-19 Rapid Test Cassette has been compared with a commercial EIA test, the results indicating high specificity and sensitivity:


IgG results:

- Relative Sensitivity: 91.8%
- Relative Specificity: 96.4%
- Accuracy: 95.0%

IgM results:

- Relative Sensitivity: 95.7%
- Relative Specificity: 97.3%
- Accuracy: 96.7%

Cat No.	Production Description	Pack Size
RADCOV1	10 Min RAPID Test Device	20 Tests
RADCOV2	10 Min RAPID Test Device with Lancet and Alcohol Swab	20 Tests
RADCOV2-02T	10 Min RAPID Test Device with Lancet and Alcohol Swab	2 Tests



Our Rapid Tests have been supplied to valued distributors from across the world to include:

ITALY | FINLAND | ROMANIA | GREECE | TURKEY | SPAIN | PORTUGAL | MEXICO | ECUADOR | PERU | COLOMBIA | BRAZIL
BULGARIA | EGYPT | ALGERIA | THE CARIBBEAN | UAE | QATAR | SAUDI ARABIA | RUSSIA | AZERBAIJAN | UNITED KINGDOM

PRESTIGE DIAGNOSTICS UK Limited

Ballymena Business Centre, Galgorm,
County Antrim, United Kingdom. BT42 1FL

www.prestigediagnosics.co.uk | info@prestigediagnosics.co.uk

T: +44 (0) 28 2564 2100



© Prestige Diagnostics U.K Ltd
August 2020 - Version 2